REF 9171/REF 9172

J-Scope™ Cannula

Tissue Visualization & Assessment Mini-Cannula Set

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
Myelotec markets this product as a sterile, single use, disposable medical device. The sterilization cycle, equipment, materials, and configuration of the device were designed for single use. Resterilization and/or re-use beyond the intent of a single use is not recommended and would pose significant risks.

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LB0308-0 2010 09

FOR USE WITH J-SCOPE™ ARTHROSCOPES

PRODUCT DESCRIPTION: Product Reference 9171 is a single use disposable device containing one (1) Mini-Cannula, one (1) Infusion Plug, one (1) Conical Obturator and one (1) Trocar. **Product Reference 9172** contains two (2) Mini-Cannulae.

INTENDED USE: When used with an endoscope, the system is designed for use as a minimally invasive, diagnostic device for visualization of joint pathology — knee, wrist, ankle, shoulder, elbow, and hips.

INDICATIONS: For introduction of fluids and instrumentation to assist in visualization of intra-articular structures for arthroscopic diagnosis, surgical and similar procedures.

- ⚠️ CONTRAINDICATIONS:** Use of the Mini-Cannula Set is contraindicated:
- for procedures involving pediatric patients.
 - when there are chronic or systematic pathologies that may increase the risk of performing arthroscopic procedures.
 - any other conditions that, in the opinion of the physician, increase the risk of performing arthroscopic surgery.

RECOMMENDED PREPARATION PRIOR TO USE:

Check the type of camera system and Arthroscope you are using. If you are in doubt about the camera system and Arthroscope compatibility, please contact your local Myelotec, Inc. agent or call Customer Service: 800.574.1633 or 770.664.4656.

If using a camera system, an appropriate video coupler with the appropriate eyepiece configuration of J-Scope Arthroscopes (either 6030-0 or 6030-30) is required.

Follow facility practice and standard procedures for ensuring proper functioning of all electrical equipment and device connectivity.

Ensure that the Arthroscope has been cleaned and sterilized per facility practice. For further information concerning the preparation, care and use of Arthroscopes, refer to the instructions that accompany the device.

Visually check the proximal window for cracks, scratches, particulate matter or sterilant build-up. Particulate matter or sterilant build-up should be carefully cleaned from the proximal window using a clean surgical swab and sterile alcohol.

⚠️ Do NOT use the arthroscope if there is any visible damage to the proximal window.

Assemble the Mini-Cannula and Arthroscope by inserting the distal tip of the arthroscope into the black scope port of the Mini-Cannula. Slowly advance the fiberoptic portion while continuously supporting the body of the arthroscope.

Lubricate the fiberoptic portion of the Arthroscope with sterile saline to ease scope insertion into the Mini-Cannula.

When fully inserted, lock the Arthroscope body in position by pushing and turning the lock pins on the body of the scope into the j-lock slots of the black scope port on the Mini-Cannula. While supporting the scope body, the yellow adjusting collar may be rotated in either a clockwise or counter-clockwise direction until the tip of the Arthroscope is flush with distal end of cannula.

⚠️ CAUTION! Forcing the arthroscope into the black scope port may kink the fiber optic arthro-scope shaft and cause irreparable damage and / or reduced optical quality.

RECOMMENDED PRESENTATION TECHNIQUE: The circulating nurse should carefully open the outer pouch and expose the contents for removal by the scrub nurse. The scrub nurse should carefully remove the inner tray from the outer tray and place it on a Mayo stand. The component tray contains the Mini-Cannula, Infusion Plug, Conical Obturator and Trocar.

Attach the light guide to the light source. Olympus , Storz , Wolf / Dyonics and ACMI light guide adapters are available from Myelotec, Inc.

Attach the light guide to the pigtail light post of the Arthroscope. Storz / Olympus and Wolf / Dyonics light post adapters are available from Myelotec, Inc. Switch on the camera and light source.

Attach the camera head to the coupler and the coupler to the Arthroscope. Connect the camera head cable to the camera control unit.

The camera should be white-balanced at this time. Ensuring adequate illumination, aim the Arthroscope at a non-reflective white background such as a surgical swab. Follow the camera manufacturer's instructions for white balancing.

Remove the Arthroscope from the Mini-Cannula by turning the body.

A small skin incision may be desirable before introducing the Mini-Cannula.

Insert the Trocar into the Mini-Cannula. The Trocar is designed to be locked into the top of the Mini-Cannula by pushing and turning the lock pins on the body of the Trocar into the j-lock slots of the black scope port. Use the completed assembly to pass through the tissue layers to the joint capsule.

Leaving the Mini-Cannula in place, remove the Trocar from the Mini-Cannula and replace it with the Conical Obturator, pushing and turning the lock pins on the body of the Conical Obturator into the j-lock slots of the black scope port. The Mini-Cannula and Conical Obturator assembly is inserted through the capsule into the joint space.

Remove the Conical Obturator leaving the Mini-Cannula in place. Insert the Arthroscope into the Mini-Cannula.

Irrigation and suction may be attached to the luer sleeve on the Mini-Cannula via the luer lock connectors. Irrigation and suction are applied alternately until a clear image is obtained on the video monitor.

If more than one Mini-Cannula is used or other portals are in use, the Infusion Plug should be used to prevent backflow of fluids.

The surgeon may now proceed with intra-articular visualization.

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